# ACCU-CHEK® Performa

## TESTS

REF 06453996 / 06454011 / 06454038 / 08966648

EN  
Device for self-testing  
Device for near-patient testing

## Intended Use

The Accu-Chek Performa system consists of the Accu-Chek Performa family of meters, Accu-Chek Performa test strips, Accu-Chek Performa control solution and Accu-Chek linearly test kit.

The test strips with the dedicated blood glucose meter are indicated to quantitatively measure glucose in fresh capillary, venous, arterial and neonatal whole blood or as aid in monitoring the effectiveness of glucose control. They are intended for self-testing by people with diabetes and for near-patient testing by healthcare professionals. They are intended for in vitro diagnostic use by healthcare professionals in clinical settings and by people with diabetes at home. Meters used in combination with an insulin pump are for home use only. For specific instructions for your meter refer to your User's Manual.

Testing sites for the Accu-Chek Performa family of meters include the finger, palm, forearm, and upper arm. Meters used in combination with an insulin pump should use fingertip testing only.

The systems are not for use in diagnosis or screening of diabetes mellitus, nor for testing neonate cord blood samples. Venous, arterial, and neonatal whole blood testing is limited to healthcare professionals only.

## Consumer Information

**Important information:** These test strips are labeled with a green  symbol to distinguish them from earlier test strips that were subject to a clinically relevant maltose interference.<sup>1</sup> The green symbol can be found on the test strip box and on the label of the test strip container.

\*Data on file

**Introduction**  
Read this package insert and the User's Manual before performing a blood glucose test.

Testing your blood glucose regularly may help you better manage your diabetes. Medical studies show that you and your healthcare professional can manage your blood glucose to near normal levels.<sup>2</sup> This can prevent or slow the development of complications from diabetes.

The package insert contains warnings and precautions:

A **WARNING** indicates a foreseeable serious hazard.  
A **PRECAUTION** describes a measure you should take to use the product safely and effectively or to prevent damage to the product.

## WARNING

**Risk of suffocation**

This product contains small parts that can be swallowed.

Keep the small parts away from small children and people who might swallow small parts.

**Risk of a serious health incident**

Failure to follow testing instructions or test strip storage and handling instructions can lead to incorrect test results.

Carefully read and follow the instructions in the User's Manual and package inserts for the test strips and control solutions.

Inspect the test strip container before using the test strips for the first time. If you see any damage to the container, such as a tear, do not use the test strips. Do not perform a control test. Contact Roche. Damaged test strips can cause inaccurate results, which could lead to improper therapy.

**Risk of infection**

Human blood is a potential source for the transmission of infection. Avoid exposing other people to contaminated components. Discard a used test strip as infectious material according to the regulations applicable in your country.

## Contents of the pack

Pack containing test strips and package inserts.

Because the reactive substances are in such small quantities, they are not considered to be hazardous materials under EU regulations. If you have any questions, contact Roche.

All components of the pack can be discarded in domestic waste.

## Test strip storage and handling

- Store the test strips at temperatures between 2 and 30 °C. Do not freeze the test strips.
- Use the test strips at temperatures between 8 and 44 °C.
- Use the test strips between 10 and 90 % humidity. Do not store the test strips in high heat and moisture areas such as the bathroom or kitchen.
- Store the unused test strips in their original test strip container with the cap closed.
- Close the test strip container tightly immediately after removing a test strip to protect the test strips from moisture.
- Use the test strip immediately after removing it from the test strip container.
- Discard the test strips if they are past the use by date. Expired test strips can produce incorrect results. The use by date is printed on the test strip box and on the label of the test strip container next to . The test strips can be used until the printed use by date when they are stored and used correctly. This applies for test strips from a new, unopened test strip container and for test strips from a test strip container already opened by the user.
- Use a test strip only once. Test strips are for single use only.

## Performing a Blood Glucose Test

Note: If your meter requires an activation chip, contact Roche to obtain one.

If you have poor circulation, testing your own blood glucose may not be right for you. Ask your healthcare professional.

Clean the puncture site before obtaining a blood drop.

- Wash your hands in warm, soapy water. Rinse and dry completely.
- Prepare the lancet device.
- Check the use by date on the test strip container. Do not use test strips past the use by date.
- Insert the test strip into the meter in the direction of the arrows. The meter turns on.
- Obtain a blood drop using the lancet device.
- Touch the blood drop to the **front edge** of the yellow window of the test strip. Remove your finger from the test strip when  appears. Do not put blood on top of the test strip.
- Remove and discard the used test strip.

Note: If the control bottle symbol and the flashing  appear on the display with your blood glucose test result, an error has occurred.

Do not act on the blood glucose result. Repeat the blood glucose test with a new test strip.

## Alternative Site Testing

You have the option of obtaining a blood sample from other sites on your body besides the fingertip. Alternative sites include the palm, forearm, and upper arm, if you are interested in alternative site testing (AST), talk to your healthcare professional first. Additional information on how to conduct AST and its limitations may be found in the User's Manual.

If you use your meter in combination with an insulin pump, only use fingertip testing.

## WARNING

**Risk of a serious health incident**

Your blood glucose level changes faster in your fingertip and palm than in the AST site (forearm and upper arm). Performing a blood glucose test with blood from the forearm or upper arm may cause you to misinterpret your actual blood glucose level, leading to improper therapy.

- Do not use AST to calibrate a continuous glucose monitoring system.
- Do not use AST to make insulin dosing calculations.
- Alternative Site Testing should be done only during steady-state times (when glucose is not changing rapidly).

## Understanding Test Results

The normal fasting glucose level for a non-diabetic adult is below 100 mg/dL (5.6 mmol/L).<sup>3</sup> The normal glucose level for a non-diabetic adult 2 hours post-meal, e.g. simulated by 75 g Oral Glucose Tolerance Test (OGTT), is less than 140 mg/dL (7.8 mmol/L).<sup>4</sup> A criterion for the diagnosis of diabetes in adults is a fasting glucose level of 126 mg/dL (higher) or, if higher, a 2-hour glucose level of 200 mg/dL (higher).<sup>5,6</sup> Results with a fasting glucose level between 100 and 125 mg/dL (5.4–6.9 mmol/L) are defined as having impaired glucose (prediabetes).<sup>5</sup> Other diagnostic criteria for diabetes exist. Consult your healthcare professional to determine if you have diabetes or not. For people with diabetes: Consult your healthcare professional for the blood glucose range appropriate for you. You should treat your low or high blood glucose as recommended by your healthcare professional.

For information on the effects and prevalence of diabetes in your area, visit the International Diabetes Federation website at [www.idf.org](http://www.idf.org) and email to [info@idf.org](mailto:idf@idf.org). For further advice or helpful information, refer to the national diabetes organization for your country.

## Unusual test results

If  is displayed on the meter, your blood glucose may be below 10 mg/dL (0.6 mmol/L).

If  is displayed on the meter, your blood glucose may be over 600 mg/dL (33.3 mmol/L).

If you receive an E-3 error message refer to your User's Manual.

## Reagent composition<sup>7</sup>

Mediator	6.72 %
Quinoprotein glucose dehydrogenase <sup>8</sup>	15.27 %
Pyruvatoquinone quinone	0.14 %
Buffer	34.66 %
Stabilizer	0.54 %
Non-reactive ingredients	42.66 %

• Minimum at time of manufacture

• From *Acacia calcicetus*, recombinant in *E. coli*, detailed description in patent application WO 2007/118647 (as "mutant 31" in table 4)

**Note:** For explanation of symbols used and a list of references, refer to the end of this package insert.

Control and linearity test kits (if available)

Accu-Chek Performa control solution — Refer to the control solution package insert for details.

Accu-Chek Performa test kit — Refer to the linearity test kit package insert for details.

Visit our website at [www.accu-chek.com](http://www.accu-chek.com) or contact the local Roche representative for more information.

## Reporting of Serious Incidents

For a patient/user/third party in the European Union and in countries with identical regulatory regime, if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and to your national authority.

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